
A Novel Medium for Increased Cytotoxicity and Serum-Minimal ex vivo Expansion of Natural Killer (NK) Cellular Immunotherapies

ABSTRACT

Ilular immunotherapy remains unacceptably variable in performan ream adoption ⁴. These hurdles are particularly true for new cancer immunotherapies w K cells ex vivo is challenging due to short half-lives, reduced functionality, and increased f serum introduces inconsistency and increases cost ⁶. Here we introduce a novel prototype mediu hich uses a unique combination of plant extracts and molecular ingredients to increas ance & proliferation, and reduce NK cell exhaustion while allowing the reduction of serum by up to 90 percer Methods: The KHYG-1 NK cell line, or primary NK cells, were brought to 2% or 2.5% serum, respectively, in the prototype media. DMEM/F12 with 20% serum (for KHYG-1) or 10% human serum (for primary cells) was the control. All media was supplemented with 100 U/mL IL-2. Cell numbers were periodically assessed, using a cell counter or flow cytometry. For cytotoxicity assays, K562 target cells were co-cultured with KHYG-1 cells (Effector:Target, E:T ratio 20:1) for 5h, or with primary NK cells (E:T = 5:1) for 18h. Cells were stained with Annexin-V and Propidium lodide to determine the levels of apoptosis and necrosis by flow cytometry. Levels of Biomarkers (e.g. CD56) were determined by Luminex (soluble biomarkers) and flow cytometry (cell-surface biomarkers). Results: Following adaptation to 2% serum, growth of KHYG-1 cells in the prototype media stayed nearly on pace with the control for over 3 weeks. Human primary NK cells grew 10-fold over 10 days in the prototype media with 2.5% serum, without feeder cells; in contrast, primary cells did not proliferate in control medium. Beta-testing with collaborating labs confirmed that the prototype media improved the proliferation of primary NK cells over a 12-day period, relative to prevalent gold standard media and methods. After culture in prototype media, both KHYG-1 cells and primary NK cells exhibited higher cytotoxic activity toward K562 and other cancer cells, compared with control medium. Secretion of interferon- γ , perforin, and granzyme A was increased in KHYG-1 cells cultured in prototype media. Additionally, NK cells are retained in a more CD56-dim phenotype, along with higher NKG2D, NKp46, CD69 and CD16: biomarker trends associated with increased cytotoxicity. Cell exhaustion biomarkers are reduced <u>Conclusions</u>: The prototype media, Enable-NK[™], increases the cytotoxic activity of NK cells against cancer cells while mitigating NK cell exhaustion. It supports proliferation of primary NK cells, even at reduced serum content, and even without feeder cells.

BACKGROUND and **INTRODUCTION**

Enable-NK[™] is a novel media formulation / media supplement for Natural Killer (NK) cells which incorporates plant-based extracts along with molecular ingredients - in a unique combination statistically optimized using Fractional-factorial Design of Experiment (DoE). Enable-NK[™] was developed under the aegis of an NIH-SBIR grant, and is at an advanced prototype testing stage.



SIGNIFICANCE


Natural Killer (NK) cellular immunotherapies which offer broader cytotoxicity, better safety and improved allogeneic potential compared to Chimeric Antigen Receptor (CAR) T cellular immunotherapies. Yet, clinical response *rates* and consistency of clinical outcomes are disappointingly low with NK cellular immunotherapies. Biomanufacturing innovations - specifically, improved culture media composition(s) represent the pathway to overcoming the entrenched hurdles of low patient response rates due to suboptimal NK cell activation, intratumoral NK cell "exhaustion", inconsistent batch-to-batch quality, and prohibitive cost. Potential benefits offered by next generation media formulations are: [a] *enhanced cytotoxicity performance* of NK cellular immunotherapies against cancer cells, resulting in higher patient response-rates in ongoing clinical trials and emergent NK cellular immunotherapies [b] *reduced batch-to-batch variability and cost*, by virtue of a serum-minimal formulation [c] *increased proliferation and viability* of bioreactor expanded NK cells (both with and without feeder


cells), resulting in improved dosage and therapy outcomes [d] *mitigated "exhaustion" of NK cells* which will improve *in vivo* persistence and patient response rates,

[e] *improved transduction* of NK cells, leading to superior CAR-NK cellular immunotherapies [f] *improved cryopreservation of NK cells* which can facilitate better logistics for effective clinical use

MATERIALS and METHODS

The co-culture studies involved the use of effector Natural Killer (NK) cells either as the KHYG-1 NK cell line or as primary NK cells [procured either from American Type Culture Collection (ATCC, Cat #PCS-800-019), or from fresh blood draws from healthy adult donors, subsequently purified from density-gradient isolated PBMC (Peripheral Blood Mononuclear Cells) using CD3+ cell depletion methods. Target cancer cells in the co-culture experiments were stained with CFSE (Carboxyfluorescein succinimidyl ester), which allowed them to be differentiated from unstained effector (NK) cells. Phycoerythrin (PE) labeled Annexin V and Propidium lodide (PI) were incorporated after the co-culture incubation period to distinguish early apoptotic (Annexin V +, PI -), late apoptotic (Annexin V +, PI +), necrotic cells (Annexin V -, PI +) and healthy cells (Annexin V -, PI -) via flow cytometry. Various incubation periods and E:T (Effector : Target) ratios were employed in the experimental strategy against a range of target cancer cells.





Authors: Marc Gillig PhD¹, Graeme Lambert¹, Matthew Forsberg PhD², Donna Sonntag¹, Jason Cahoon¹, Rafet Amoor¹, Borom Chean¹, Lea Picard PhD³, Sabine Wingert³, Caroline Hull PhD⁴, Jacki Kornbluth PhD⁵, Carsten Watzl PhD³, Christian Capitini MD², Rachit Ohri PhD^{1*}. (*Corresponding Author, Email: <u>rachit@enablelifesciences.com</u>)

Affiliations: 1: Enable Life Sciences LLC, Worcester MA. 2: Univ of Wisconsin, Madison. 3: Leibniz Research Centre, Germany. 4: Leucid Bio Ltd, London, UK. 5: St. Louis University.

RESULTS organized by Performance Attribute of NK Cells

