

DEPARTMENT OF Pediatrics UNIVERSITY OF WISCONSIN School of Medicine and Public Health



LEIBNIZ RESEARCH CENTRE FOR WORKING ENVIRONMENT


ABSTRACT

cells have emerged as a promising cancer immunotherapy due to their innate abilities (e.g. antige ns) to recognize and kill tumor cells, as well as their allogeneic potential. However, NK cells are ex vivo; therefore, new cell culture media technologies are needed for ex vivo etaining or enhancing their therapeutic functions. In this study, we investigated expans and cytotoxic activity of primary human NK cell cultures using a novel 2-part media formula cular ingredients and plant-derived extracts in an optimal media recipe. We achieved 10-fc expansion in the absence of feeder cells over 14 days, with 10% FBS in the presence of IL-2 and IL-15 – much higher compared to the primary NK cells from the same donor cultured in control media conditions i.e. RPMI with 10% FBS in the presence of IL-2 and IL-15. Cells cultured in the novel media had lower expression of the inhibitory receptor NKG2A, and exhaustion markers LAG-3 and TIM-3. These cells also exhibited higher expression of CD56 after 14 days of expansion compared to the freshly isolated cells in control media, but lower expression of CD56 after 14 days of expansion in control media. Using the Incucyte real-time quantitative live-cell imaging platform, enhanced cytotoxicity was observed in novel media-expanded NK cells compared to freshly isolated NK cells against the human melanoma cell line M21 and the human pediatric neuroblastoma cell line CHLA-20. When the anti-GD2 monoclonal antibody Ch 14.18 was added to the coculture antibody dependent cellular cytotoxicity (ADCC) against the target cell lines M21 and CHLA-20 was also seen to be enhanced in novel media-expanded NK cells compared to freshly isolated NK cells from the same donor in control media. When the novel media were incorporated into feeder cell-based NK cell expansion protocols (e.g. using engineered K562 cells), the cultures using the novel media outperformed the feeder-based culture using control media in terms of both proliferation and cytotoxic activity. Finally, the cytotoxic activity of cryopreserved NK cells expanded in the novel media was superior compared to control media expanded NK cells. In summary, this novel 2-part culture medium performs well in expanding NK cells, both with and without feeder cells, and increases cytotoxic performance; the underlying biological mechanisms are manifested in a consistent expression pattern of cell surface biomarkers relevant to NK cell biology.

BACKGROUND and INTRODUCTION

Enable-NK[™] is a novel media formulation / media supplement for Natural Killer (NK) cells which incorporates plant-based extracts along with molecular ingredients - in a unique combination statistically optimized using Fractional-factorial Design of Experiment (DoE). Enable-NK[™] was developed under the aegis of an NIH-SBIR grant, and is currently available as a Research Use Only (RUO) product.

Two-Phase Media Format




SIGNIFICANCE


Natural Killer (NK) cellular immunotherapies which offer broader cytotoxicity, better safety and improved allogeneic potential compared to Chimeric Antigen Receptor (CAR) T cellular immunotherapies. Yet, clinical response *rates* and consistency of clinical outcomes are disappointingly low with NK cellular immunotherapies. Biomanufacturing innovations - specifically, improved culture media composition(s) represent the pathway to overcoming the entrenched hurdles of low patient response rates due to suboptimal NK cell activation, intratumoral NK cell "exhaustion", inconsistent batch-to-batch quality, and prohibitive cost. Potential benefits offered by next generation media formulations are: [a] enhanced cytotoxicity performance of NK cellular immunotherapies against cancer cells, resulting in higher patient response-rates in ongoing clinical trials and emergent NK cellular immunotherapies [b] reduced batch-to-batch variability and cost, by virtue of a serum-minimal formulation [c] *increased proliferation and viability* of bioreactor expanded NK cells (both with and without feeder cells), resulting in improved dosage and therapy outcomes


[d] *mitigated "exhaustion" of NK cells* which will improve *in vivo* persistence and patient response rates, [e] *improved transduction* of NK cells, leading to superior CAR-NK cellular immunotherapies [f] *improved cryopreservation of NK cells* which can facilitate better logistics for effective clinical use

MATERIALS and METHODS

The control media used throughout the experiments was either RPMI-1640 supplemented with 10%FBS, 1% L-glutamine, 1% penicillin/streptomycin or RPMI-1640 supplemented with 10% human AB serum, 1% penicillin/streptomycin, 1% non-essential amino acids, 1% HEPES buffer. Additional media was added to the cultures every 2-3 days, to keep the cell concentration between 1x10⁶/ml and 2x10^{^6}/ml. 100 IU/ml IL-2 and 10 ng/ml IL-15 was also added to the cultures every 2-3 days. The feeder cells used were K562 cells engineered to express 41BBL and IL-15 (K562-41BBL-IL15). For all Incucyte assays, 10,000 target cells (either CHLA-20-GFP or M21-GFP) were plated in each well of a 96-well plate. The media for the Incucyte assays was DMEM supplemented with 10% FBS. 1ug of the anti-GD2 monoclonal antibody Ch14.18 was added to each ADCC well. Flow cytometry was performed on day 14 or day 18 using the Attune NxT flow cytometer, with 200,000 events collected for each sample. Antibodies were purchased from Biolegend (CD3-BV605, CD56-FITC or CD56-PE-Dazzle 594, CD16-AF700, LAG-3-BV711, NKG2A-AF700, TIM-3-BV510).







Novel NK cell culture medium supports enhanced proliferation and cytotoxic activity of peripheral blood NK cells.

<u>Authors</u>: FORSBERG Matthew¹, SCHUERMAN Elizabeth¹, PICARD Lea², WINGERT Sabine², GILLIG Marc³, OHRI Rachit ^{3*}, WATZL Carsten², CAPITINI Christian¹

(*Corresponding Author, Email: <u>rachit@enablelifesciences.com</u>)

Affiliations: 1: Univ of Wisconsin, Madison. 2: Leibniz Research Centre, Germany. 3: Enable Life Sciences LLC, Worcester MA.

RESULTS organized by Performance Attribute of NK Cells

